


# How to Clean & Store Your Protective Garments

Maintain your investment and safeguard your health

Improper storage and cleaning can degrade and compromise your protective garments. Here are some tips to get the maximum protection for as long as possible.



## CLEANING

### Do's

- ✓ Spot clean with a damp cloth and dishwashing detergent. Dry with a cloth.
- ✓ If you need to scrub, only use very soft bristles like a toothbrush.
- ✓ Periodically sanitize with a QACs cleaning agent such as DisCide.

### Do not's

- ✗ Do not soak or submerge in water.
- ✗ Do not use washers or dryers.
- ✗ Do not autoclave.
- ✗ Avoid harsh cleaning agents like bleach.

## STORAGE

### Do's

- ✓ Hang aprons and vests by the shoulders or straps using an apron hanger.
- ✓ While transporting, roll and place in an apron bag with a tube.

### Do not's

- ✗ Avoid folding or creasing the garment.
- ✗ Do not hang by a hook.
- ✗ Avoid extreme heat, cold or humidity.



Protective garments are an investment in your health.  
Take care of them and they will take care of you!