

Fully Customized and Made to Order Table Pads and Covers

Increase patient safety and reduce length of stay with high-quality, extremely durable, and economically advantageous table pads and covers from Burlington Medical.



Dedicated account managers who provide the highest level of customer service and product knowledge.



Top quality products and more than 40 years of manufacturing experience.



Fastest order turn-around in the industry.

TO PLACE YOUR ORDER

Talk to your representative | Call 800-221-3466

BurMed.com | 1-800-221-3466 | 3 Elmhurst St. / PO BOX 3194, Newport News, VA 23603



COATED RADIOLOGY PADS

Burlington Medical table pads are available in any configuration—we can fit any table pad for radiology, cath lab, OR, and stretcher and gurney pads.

All table pads come with a two-year manufacturer's warranty.



CATH LAB | IR | EP

Pick your own configuration:

Rectangle | 2-Step | 3-Step (pictured)

Tri Foam is the recommended foam choice in these pads for superior pressure management during long cases

OPERATING ROOM

Fit table pads of any size, shape, manufacturer, and model:

3-Piece (pictured) | 4-Piece | Torso Cut-out

Sectional replacements are available, so you buy only the pieces you need

STRETCHER & GURNEY

- Customizable to fit all popular stretcher and gurney models
- Can have angled corners or large rounded corners
- Optional 18" hook-and-loop strips
- Available in a range of foam and covering configurations to meet your specifications

RADIOLOGY

COVERED RADIOLOGY TABLE PADS

- A range of configuration options to meet the needs of your specific application

COATED RADIOLOGY TABLE PADS

- Stain resistant
- Extremely durable
- Standard full length size: 24" x 77" x 1"



CUSTOM-MADE TABLE PADS

1 CUSTOM MEASUREMENTS

Your Account Manager will measure your pad/table for proper dimensions.

2 FOAM CONFIGURATION

Ester Foam

Economical open-cell foam that offers adequate cushioning for healthy adults over short periods of time. Radiolucent for appropriate procedures.

Depths: 1", 2", 3", 4", or custom

C Foam

Dense closed-cell foam that yields to the touch but won't bottom out. Radiolucent for appropriate procedures.

Depths: 0.5", 1", 1.5", or custom

Q Foam

Memory Foam and C Foam Bi-layer follow the contour of the body, reducing pressure points. Good for high-risk patients and long procedures. Radiolucent for appropriate procedures.

Depths: 1.5", 2", or custom

Tri Foam

Combines the benefits of C Foam, Ester, and Memory Foam to create a maximum comfort table pad. Considerably reduces pressure points and is suitable for long procedures. Radiolucent for appropriate procedures.

Depths: 2.5", 3", or custom

3 COVER MATERIAL AND COLOR

Choose the perfect covering for your table pad to increase comfort, resist fluid, or cut down on static electricity.

Dartex™ has been the world leader in technical support surfaces for pressure injury prevention for over 40 years. These fabrics offer unmatched durability, fluid resistance, and conformability to reduce pressure on patients' skin surfaces and offer a simple wipeable surface that can stand up to 10,000 ppm bleach when applied for a limited time. Anti-fungal properties also protect the effectiveness of the coating for the lifetime of the fabric.

Available in the following colors

#801
Black

#802
Navy

#803
Green

#804
Orange

#805
Red

#806
Purple

#807
Pink

#808
Blue

Dartex Conductive #850 has the same great benefits as standard Dartex fabrics with built-in anti-static technology that minimizes electrical interference with sensitive equipment.

Available in the following color

#850
Black

Vinyl is a heavyweight, durable, non-elastic vinyl commonly used on table pads. Best for pads that are exposed to a lot of fluid. Not recommended with Q Foam or Tri Foam.

Available in the following colors

#302
Royal

#306
Black

#996
Grey

4 DESIRED FEATURES

Mega Grip is a performance fabric for the bottom surface of your table pads that boasts an impressive coefficient of friction on almost all surfaces. Avoid the hassles and infection control concerns of hook-and-loop and upgrade to a better solution for your table pads.

Closures determine how you open the cover to access the table pad, or sew it shut.

Available in the following

Zipper | Hook-and-Loop | Sewn Shut



Outlast® is a high performance phase-change materials improve the temperature regulation on a variety of applications. Create a more comfortable and stable temperature environment on your table pads.

Other available features

Hanging Tabs | Handles | Restraints