

Made to Order Table Pads and Covers

Increase patient safety and reduce length of stay with high-quality, extremely durable, and economically advantageous table pads and covers.



Fastest turn-around in the industry



Decades of expertise for improved patient outcomes



Our dedicated Account Managers provide the highest level of customer service and product knowledge.



1 CUSTOM MEASUREMENTS

Your Account Manager will measure your pad/table for proper dimensions. In addition to custom solutions, we also offer a wide selection of standard pads designed to fit many makes and models, including the following table brands: AMSCO, Steris, Sktryer, Skytron, Hilrom and Shampaine.



SCAN OR CLICK TO
VIEW STANDARD
TABLE PAD OPTIONS

2 FOAM CONFIGURATION

All are radiolucent for appropriate procedures.

Ester Foam

Economical open-cell foam that offers adequate cushioning for healthy adults over short periods of time.

Depths: 1", 2", 3", 4"

C Foam

Dense closed-cell foam that yields to the touch but won't bottom out.

Depths: 0.5", 1", 1.5"

Q Foam

Memory Foam and C Foam Bi-layer follow the contours of the body, reducing pressure points. Good for high-risk patients and long procedures.

Depths: 1.5", 2"

Tri Foam

Combines the benefits of C Foam, Ester, and Memory Foam to create a maximum comfort table pad. Considerably reduces pressure points and is suitable for long procedures.

Depths: 2.5", 3"

4 OPTIONAL FEATURES

Mega Grip is a performance fabric for the bottom surface of your table pads that boasts an impressive coefficient of friction on almost all surfaces. Avoid the hassles and infection control concerns of hook-and-loop and upgrade to Mega Grip.

Closures determine how or if the table pad is accessible.

Available in the following styles:

Zipper | Hook-and-Loop | Sewn Shut

3 COVER MATERIAL AND COLOR

Choose the perfect covering for your table pad to increase comfort, resist fluid, or cut down on static electricity.

Dartex™ fabrics offer unmatched durability, fluid resistance, anti-fungal properties, and conformability. They reduce pressure on patients' skin surfaces with a simple wipeable surface that can stand up to 10,000 ppm bleach when applied for a limited time.

Available in the following colors:

#570 Dartex Performance Black	#571 Dartex Performance Jade	#572 Dartex Performance Red	#573 Dartex Performance Purple
#574 Dartex Performance Blue	#575 Dartex Performance Navy	#803 Dartex Care Green	#804 Dartex Care Orange

Dartex Conductive #850 has the same great benefits as standard Dartex fabrics with built-in anti-static technology that minimizes electrical interference with sensitive equipment.

Available in the following color:

#850
Black

Vinyl is heavyweight, durable, and non-elastic. Best for pads that are exposed to a lot of fluid.

Not recommended with Q Foam or Tri Foam.

Available in the following colors:

#302 Royal	#306 Black	#996 Grey
----------------------	----------------------	---------------------



Outlast® is a high performance phase-change material that improves the temperature regulation on a variety of applications. Create a more comfortable and stable temperature environment on your table pads.

Other Available Features:

Hanging Tabs | Handles | Restraints